

ReSound Vida

Rechargeable Questions and Answers



<p>What is included in the ZPower kit?</p>	<p>The ZPower kit contains two battery door modules, two silver-zinc rechargeable batteries and a charger.</p>
<p>Is there a ZPower kit for monaural users?</p>	<p>No, ZPower does not have a kit for monaural users. A monoraul fitting will need to have an order for a battery door, silver-zinc rechargeable battery and a charger separately.</p>
<p>What is the life of a silver-zinc rechargeable battery?</p>	<p>ZPower silver-zinc rechargeable batteries are designed to last more than a year based on daily wear.</p>
<p>What if members use a silver-zinc rechargeable battery longer than one year?</p>	<p>Members will still be able to use their ReSound Vida Rechargeable hearing aids as normal, they may just notice the rechargeable battery life becoming progressively shorter.</p>
<p>How many hours of use can members expect from a single charge?</p>	<p>ReSound Vida Rechargeable hearing aids deliver continuous, uninterrupted power up to 16 hours with no streaming. They will last up to 14 hours with 2 hours of direct streaming.</p>
<p>How long does it take to charge the battery?</p>	<p>The battery will be fully charged in seven hours or less. Charging the hearing aids every day is recommended.</p>

What if a ZPower silver-zinc battery is used without a ZPower battery door module?	A ZPower battery door module is required to use ZPower silver-zinc batteries. Without ZPower battery door modules, ZPower silver-zinc batteries will not power the hearing aids.
Can disposable zinc-air batteries be used in a ReSound Vida Rechargeable hearing aid?	Yes, size 312 zinc-air batteries are compatible.
What happens if members accidentally attempt to charge a disposable zinc-air battery in the charger?	The ZPower charger is able to distinguish between rechargeable silver-zinc batteries and disposable zinc-air batteries. A red LED light on the charger will illuminate if members attempt to charge disposable batteries.
Are the ZPower battery door modules and silver-zinc rechargeable batteries nanocoated?	ZPower battery door modules are nanocoated (not P2i) and ZPower silver-zinc batteries are not. However, silver-zinc batteries are sealed so if a member drops them in water, the batteries can still be used once completely dried.
Are there any changes to the warranty once a ReSound Vida Rechargeable hearing aid is modified with a ZPower battery door?	No, the product warranty remains the same. The ZPower charger and battery door modules are covered on a separate accessory warranty.

Manufacturer according to FDA:

GN ReSound North America
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-882-3636
future-resound.com

Manufacturer according to Health Canada:

ReSound Canada
303 Supertest Road
Toronto, ON M3J 2M4
Canada
1-888-737-6863
future-resound.com



MK605078 Rev A 2018.06

